

Supplemental Material. Study protocol. The Dexcom G6-Pro accuracy was assessed over the 10 days of sensor wear by comparing time matched values to:

- i) Self-monitoring of blood glucose (SMBG) performed by the participants 4 to 7 times per day while outside the hemodialysis sessions (blue boxes).
- ii) Venous blood glucose measurements were obtained during hemodialysis sessions on Days 4, 6 and 8 (orange boxes). Blood samples from the venous line of the hemodialysis circuit were obtained by study staff at pre-specified time points to capture the early shifts in glucose during the hemodialysis session.

